

# Certificate of Analysis

No. CA22229

Product: Eichrom's LN-2 Resin Cartridges  
Lot #: L2SR22B  
Date of Manufacture: 7-OCT-2022

Eichrom's LN-2 Resin cartridge Lot L2SR22B has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, # L2S100322, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

| <u>Test Value</u> | <u>Specification</u> | <u>Result</u> |
|-------------------|----------------------|---------------|
| Appearance        | White Beads          | Pass          |
| Lu Peak Maximum   | 11 - 15 BV           | 11.9 BV       |
| Yb Peak Maximum   | 6.0 - 9.0 BV         | 7.0 BV        |
| Seperation Factor | $\geq 1.6$           | 1.8           |

  
Jill Bryant  
Quality System Coordinator

October 10, 2022