

### Certificate of Analysis


No. CA22235

Product: Eichrom's UTEVA® Resin 2mL Columns  
Lot #: UT22J  
Date of Manufacture: 11-OCT-2022

Eichrom's UTEVA® Resin 2mL Columns Lot UT22J has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, UTA-MAE-189-89, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

### ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.1 N nitric acid	>100	167.4
% <sup>230</sup> Th stripped in 5M HCl/0.05M oxalic acid		
Least acceptable mean (4 rep)	90.0	102.2
Max. standard deviation	10.0	3.0
% <sup>230</sup> Th stripped in 1M HCl		
Maximum acceptable mean (5 rep)	0.5	0.1
Max. standard deviation	1.0	0.0
% <sup>233</sup> U stripped in 1M HCl		
Least acceptable mean (5 rep)	90.0	101.4
Max. standard deviation	10.0	4.8
% <sup>233</sup> U stripped in 5M HCl/0.05M oxalic acid		
Maximum acceptable mean (5 rep)	1.0	0.2
Max. standard deviation	1.0	0.0

  
Sarah McAlister  
Quality Manager

13 Oct 2022  
October 13, 2022