

# Certificate of Analysis

No. CA 22292

Product: Eichrom's Anion Exchange Resin QML Cartridges

Lot #: A8MR22D


Date of Manufacture: 14-DEC-2022

Eichrom's Anion Exchange resin QML cartridge Lot A8MR22D has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, 070522A8M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	0.03
Wet Volume Capacity	≥ 1.2 meq/mL	1.4

Jill Bryant  
Quality System Coordinator



December 19, 2022