

Certificate of Analysis

No. CA 23019

Product: SCX Cartridges, 1mL
Part No: SCX1ML-R50
Lot No: SCXR23A
Manufactured: 26-JAN-2023
Date:

Eichrom's SCX 1mL Cartridge Lot # SCXR23A has been manufactured according to Eichrom's standard operating procedures. The silica used to manufacture this lot, SCX060922, has been tested in accordance with Eichrom's quality assurance procedures and purity testing. This lot of product meets all Eichrom quality specifications.

No products of animal or human origin are used in the manufacture of this product.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Recovery of Ga (%)	>98.0	101.0

Jill Bryant
Quality System Coordinator



January 31, 2023