

**Certificate of Analysis**  
No. CA23125

Product: Eichrom's DGA-Branched Resin  
Lot #: DBS060623  
Date of Manufacture: 6-JUN-2023

Eichrom's DGA-Branched Resin Lot # DBS060623 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% <sup>241</sup> Am stripped in 0.01M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	90.0	98.2
Max. standard deviation	2.5	0.2
% <sup>241</sup> Am stripped in 0.5M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	0.0
Max. standard deviation	1.0	0.2
% <sup>233</sup> U stripped in 0.5M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	95.0	98.3
Max. standard deviation	2.5	0.2
% <sup>233</sup> U stripped in 0.01M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	1.0
Max. standard deviation	1.0	0.2

Jill Bryant  
Quality System Coordinator

June 14, 2023