

Certificate of Analysis


No. CA23137

Product: Eichrom's LN-2 Resin
Lot #: L2S061923
Date of Manufacture: 19-JUN-2023

Eichrom's LN-2 Resin Lot # L2S061923 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
Lu Peak Maximum	11 - 15 BV	12.8 BV
Yb Peak Maximum	6.0 - 9.0 BV	7.9 BV
Seperation Factor	≥ 1.6	1.7



Jill Bryant
Quality System Coordinator

July 6, 2023