

Certificate of Analysis

No. CA23163

Product: Eichrom's LN Resin
Lot #: LNF072523
Date of Manufacture: 26-JUL-2023

Eichrom's LN Resin Lot # LNF072523 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2760



Jose A Munoz
Quality Analyst

July 31, 2023