

# Certificate of Analysis

No. CA 23087

Product: Eichrom's Nickel Resin  
Lot #: NIA050223  
Date of Manufacture: 2-MAY-2023

Eichrom's Nickel Resin Lot # NIA050223 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Ave. % Ni Recovery	80-105%	101.8 %
Ave. % Ni in Load/Rinse	<1%	0.2 %
Ave. % Ni in Strip	>99%	99.8 %

  
Jill Bryant  
Quality System Coordinator

May 8, 2023