

Certificate of Analysis

No. CA 23302

Product: Eichrom's Anion Exchange Resin Cartridges

Lot #: A8MR23C

Date of Manufacture: 15-DEC-2023

Eichrom's Anion Exchange resin cartridge Lot A8MR23C has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, 063021A8M, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White to light Yellow	Pass
Extractables	< 1 mg (residue) / g (resin)	n.d.
Wet Volume Capacity	≥ 1.2 meq/mL	1.3


Jill Bryant
Quality System Coordinator

December 18, 2023