

Certificate of Analysis

No. CA 23278

Product: Eichrom's Actinide Resin Cartridges
Lot #: ACSR23E
Date of Manufacture: 6-NOV-2023

Eichrom's Actinide Resin cartridge Lot # ACSR23E has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, ACS080923, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.5 M HCl acid	≥ 300	557



Jose Munoz
Quality Analyst

November 8, 2023