

**Certificate of Analysis**  
No. CA23232

Product: Eichrom's DGA-Branched Resin Cartridges  
Lot #: DBSR23F  
Date of Manufacture: 28-SEP-2023

Eichrom's DGA-Branched Resin cartridge Lot DBSR23F has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, DBS090523, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
% <sup>241</sup> Am stripped in 0.01M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	90.0	97.0
Max. standard deviation	2.5	0.2
% <sup>241</sup> Am stripped in 0.5M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	0.6
Max. standard deviation	1.0	0.3
% <sup>233</sup> U stripped in 0.5M HNO <sub>3</sub>		
Least acceptable mean (4 rep)	95.0	97.7
Max. standard deviation	2.5	0.1
% <sup>233</sup> U stripped in 0.01M HNO <sub>3</sub>		
Maximum acceptable mean (4 rep)	2.0	1.61
Max. standard deviation	1.0	0.1



Jose Munoz  
Quality Analyst

September 29, 2023