

Certificate of Analysis

No. CA23172

Product: Eichrom's DGA-Normal Resin HML Cartridges
 Lot #: DNSR23F
 Date of Manufacture: 2-AUG-2023

Eichrom's DGA-Normal Resin HML Cartridge Lot DNSR23F has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, DNS052123, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution ratio, Eu in 0.1M nitric acid	≥900	1951
% ²⁴¹ Am stripped in 0.1M HCl		
Least acceptable mean (4 rep)	97.0	99.0
Max. standard deviation	2.5	0.3
% ²⁴¹ Am stripped in 0.1M HNO ₃		
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.1
% ²³³ U stripped in 0.1M HNO ₃		
Least acceptable mean (4 rep)	97.0	98.5
Max. standard deviation	2.5	0.1
% ²³³ U stripped in 0.1M HCl		
Maximum acceptable mean (4 rep)	1.0	0.7
Max. standard deviation	0.5	0.1



Sarah McAlister
Quality Manager

August 2, 2023