



Certificate of Analysis

No. CA23223

Product: Eichrom's LN-3 Resin Cartridges
Lot #: L3SR23A
Date of Manufacture: 18-SEP-2023

Eichrom's LN-3 Resin cartridge Lot L3SR23A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, L3S030622, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
La in Load/Rinse (%)	95-105	99.9
Be in Load/Rinse (%)	95-105	99.5
U in Load/Rinse (%)	< 3	0.2
La in HEDPA Strip (%)	< 1	0.1
Be in HEDPA Strip (%)	< 1	0.5
U in HEDPA Strip (%)	95-105	99.8


Jill Bryant
Quality System Coordinator

September 19, 2023