

# Certificate of Analysis

No. CA23267

Product: Eichrom's LN-3 Resin Cartridges  
Lot #: L3SR23B  
Date of Manufacture: 25-OCT-2023

Eichrom's LN-3 Resin cartridge Lot L3SR23B has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, L3S030622, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
La in Load/Rinse (%)	95-105	99.9
Be in Load/Rinse (%)	95-105	99.5
U in Load/Rinse (%)	< 3	0.2
La in HEDPA Strip (%)	< 1	0.1
Be in HEDPA Strip (%)	< 1	0.5
U in HEDPA Strip (%)	95-105	99.8

  
Jill Bryant  
Quality System Coordinator

October 26, 2023