

# Certificate of Analysis

No. CA23207


Product: Eichrom's LN Resin Columns  
Lot #: LN23A  
Date of Manufacture: 30-AUG-2023

Eichrom's LN Resin column Lot LN23A has been manufactured according to Eichrom standard operating procedures. The resin used to manufacture this lot, LNA051822, has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2086

Jill Bryant  
Quality System Coordinator



August 31, 2023