

Certificate of Analysis

No. CA23194

Product: Eichrom's LN Resin
Lot #: LNA082423
Date of Manufacture: 24-AUG-2023

Eichrom's LN Resin Lot # LNA082423 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2621



Jill Bryant
Quality System Coordinator

August 24, 2023