



# Certificate of Analysis

No. CA23269

Product: Eichrom's LN Resin  
Lot #: LNF103123  
Date of Manufacture: 31-OCT-2023

Eichrom's LN Resin Lot # LNF103123 has been manufactured according to Eichrom standard operating procedures. This lot has been tested in accordance with Eichrom quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

## ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 0.1 N nitric acid	>700	2866.5

Jose Munoz  
Quality Analyst

October 31, 2023