

**Certificate of Analysis**  
No. CA 24026

Product: Eichrom's TEVA Resin  
Lot #: TES020724  
Date of Manufacture: 7-FEB-2024

Eichrom's TEVA Resin Lot # TES020724 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	400
% <sup>230</sup> Th stripped in 6M HCl		
Least acceptable mean (4 rep)	93.0	98.6
Max. standard deviation	2.5	0.1
% <sup>230</sup> Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.1
Max. standard deviation	1.5	0.1
% <sup>239</sup> Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	99.4
Max. standard deviation	2.5	0.6
% <sup>239</sup> Pu stripped in 6M HCl		
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.2

  
Jose Munoz  
Quality Analyst

February 13, 2024