



eichrom
TECHNOLOGIES

A GCI COMPANY

Certificate of Analysis
No. CA24043

Product: Eichrom's TEVA Resin 1mL Cartridges
Lot #: TESR24B
Date of Manufacture: 26-Feb-2024

Eichrom's TEVA Resin 1mL Cartridges Lot TESR24B has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, TES072723, has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

ANALYSIS RESULTS

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	428.2
% ²³⁰ Th stripped in 6M HCl		
Least acceptable mean (4 rep)	93.0	98.9
Max. standard deviation	2.5	0.1
% ²³⁰ Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.1
Max. standard deviation	1.5	0.1
% ²³⁹ Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	96.8
Max. standard deviation	2.5	1.5
% ²³⁹ Pu stripped in 6M HCl		
Maximum acceptable mean (4 rep)	1.0	0.1
Max. standard deviation	0.5	0.0

Jill Bryant
Quality System Coordinator

February 28, 2024

Eichrom Technologies LLC
1955 University Lane, Lisle, IL 60532, USA
www.eichrom.com

npo

eichrom