

Certificate of Analysis

No.

Product:

Part Number:

Lot Number:

Manufactured Date:

This product has been manufactured according to Eichrom's standard operating procedures. The resin used to manufacture this lot, , has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

Analysis

<u>Test</u>	<u>Specification</u>	<u>Result</u>
Appearance	White Beads	Pass
Distribution Ratio, Eu in 0.1M HNO ₃	≥ 900	
% ²³³ U stripped in 0.1M HCl		
Max. acceptable mean	1.0	
Max. standard deviation	0.5	
% ²⁴¹ Am stripped in 0.1M HCl		
Min. acceptable mean	97.0	
Max. standard deviation	2.5	
% ²³³ U stripped in 0.1M HNO ₃		
Min. acceptable mean	97.0	
Max. standard deviation	2.5	
% ²⁴¹ Am stripped in 0.1M HNO ₃		
Max. acceptable mean	1.0	
Max. standard deviation	0.5	