

**Certificate of Analysis**  
No. CA 24152

Product: Eichrom's TEVA Resin  
 Lot #: TEA071124  
 Date of Manufacture: 11-JUL-2024

Eichrom's TEVA Resin Lot # TEA071124 has been manufactured according to Eichrom's standard operating procedures. This lot has been tested in accordance with Eichrom's quality assurance policy and procedures. This lot of product meets all Eichrom quality specifications.

**ANALYSIS RESULTS**

<u>Test Value</u>	<u>Specification</u>	<u>Result</u>
Appearance	White beads	Pass
Distribution Ratio, 3.0 N nitric acid	>200	468.6
% <sup>230</sup> Th stripped in 6M HCl		
Least acceptable mean (4 rep)	93.0	97.8
Max. standard deviation	2.5	0.2
% <sup>230</sup> Th stripped in 0.02M HCl + 0.02M HF		
Maximum acceptable mean (4 rep)	2.0	0.7
Max. standard deviation	1.5	0.2
% <sup>239</sup> Pu stripped in 0.02M HCl + 0.02M HF		
Least acceptable mean (4 rep)	90.0	97.38
Max. standard deviation	2.5	1.6
% <sup>239</sup> Pu stripped in 6M HCl		
Maximum acceptable mean (4 rep)	1.0	0.0
Max. standard deviation	0.5	0.1

  
 Jose Muñoz  
 Quality Analyst

July 16, 2024