

TYPE A PACKAGE FOR 11 SPECT UNIT DOSE SYRINGES

Model: RSS SPECT - UNIT DOSE SHIELD (UDS)

Eichrom Technologies designs and manufactures compact Type A shipping systems for the safe transport of radionuclides. Our shielded RAD SHIPPING SYSTEMS are built for strength, efficiency, and ergonomic use, ensuring secure transport, easy handling, and full regulatory compliance. We offer both standard and customizable solutions with fast delivery and expert support. To place an order, request pricing, or ask questions, contact us today. ☎ +1 (630) 963-0320 ✉ info@eichrom.com

- ✓ **1/4-turn syringe pig (UDS) twist-off** reduces handling time and radiation exposure
- ✓ **Universal design** allows new products to integrate with the current inventory
- ✓ **Thick, heavy-duty canvas case** is flexible, lightweight, strong, and durable
- ✓ **Hygienic UDS units** are easy to clean and dishwasher-safe



Durable stainless steel syringe pigs (UDS)
won't crack, break, or open when dropped

Universal design
allows new products
to integrate with the
current inventory



Sealed O-ring design eliminates
the need for absorbent material

Triple-start threads for efficient 1/4-turn operation enable quick opening
and closing to minimize radiation exposure and reduce handling time

Perfectly flat bottom ensures stability when upright to
prevent spilling radioactive material or product loss



Removable clear envelopes
attach securely with Velcro

Thick, heavy-duty canvas case
is lightweight, strong, and durable, with reinforced corners and a sturdy bottom



Durable construction features a smooth-glide nylon zipper that won't rust or corrode, reinforced rivets, and strong woven handles

RSS SPECT – UNIT DOSE SHIELD SPECIFICATIONS

Part Number	SH-SYS-UDS-SPECT
Description	Rad Shipping System for (11) Spect Unit Dose Shields, Includes Spect Case and (11) Spect Unit Dose Shields
Authorized Contents	(11) 3 - 10mL Syringes
Case Color	Black
Total Weight	42 lb (19.1 kg)
Max Activity F-18 for Yellow II (Approximate)	N/A

REGULATORY COMPLIANCE



DOT Type A Packaging
Requirements



IATA Dangerous Goods
Regulations



IAEA Regulations For The Safe
Transport Of Radioactive Material

MORE INFO?

Contact us today:
+1 (630) 963-0320
info@eichrom.com

